



Salient Features

- 1. Fully automatic untouched process for liquid oral manufacturing.
- 2. cGMP compliance with CIP/SIP facility.
- 3. Capacities available from 100 ltrs. to 20000 ltrs. batch capacity
- 4. Fully automatic plant consisting of Sugar syrup vessel, syrup manufacturing vessel and storage vessel with interconnecting pipe lines, transfer pumps and integrated control panel.
- 5. Working platform for 1000 ltrs. and above capacity plants.
- 6. Entire plant can be operated through integrated control panel, by 01 skilled & 01 semi skilled operator.
- 7. Sugar charging by means of vacuum transfer or mechanical charging can be provided optionally.
- 8. All interconnecting pipe lines are internally electro polished and with easy openabletriclover fittings.
- 9. All vessels are equipped with bottom entry agitator, which leads to less maintenance, lower energy consumption and better mixing of products.
- 10. Different plant layouts to suit site conditions or customers' requirement.
- II. Cleaning validation by SWEB testing match international standard less than 10 PPM OR by CIP system with online conductivity.
- 12. Load Cell for accurate weighing of final products.



Salient Features

- 1. Fully automatic untouched process for Oingment / Cream / Lotion manufacturing.
- 2. cGMP compliance with CIP/SIP facility.
- 3. Capacities available from 25 kgs. to 2500 kgs. batch capacity
- 4. Fully automatic plant consisting of Wax Phase vessel, Water phase vessel, manufacturing vessel & storage vessel with interconnecting pipe lines, transfer pumps and integrated control panel.
- 5. Working platform for 200 kgs. and above capacity plants.
- 6. Entire plant can be operated through integrated control panel, by 01 skilled & 01 semi skilled operator.
- 7. Wax & water phase transferred by means of vacuum transfer, to main manufacturing vessel.
- 8. Semi contra type mixer for cream / ointment manufacturing. Contra rotary Mixer or Planetary Mixer can be provided upon demand.
- 9. Main mixer shall be provided with heating / cooling jacket with digital temperature controller.
- 10. Plants available for sterile ointments manufacturing.
- 11. Process piping is highly Electro polished and joints are DIN/TC standard for easy opening & re-fixing.
- 12. Hydraulic Power Pack with cylinder to lift complete agitator assembly with top lid for product recovery / Easy to cleaning / Change over of batch & for maintenance.
- 13. VFD for Anchor agitators to vary the RPM as per product / Recipe.
- 14. Different plant layouts to suit site conditions or customers' requirement.
- 15. Cleaning validation by SWEB testing match international standard less than 10 PPM OR by CIP system with online conductivity.
- 16. Load Cell for accurate weighing of final products.
- 17. Vacuum system save timing to transfer product from one vessel to another vessel and zero retention in pre phase vessel & inter connecting pipeline.

